Instructions For Use

REF: PMP 802 / PMP 804 / PMP 803

Chromoprobe Multiprobe® OctoChrome™

FOR PROFESSIONAL USE ONLY

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Further information available at www.cytocell.com

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei from fixed cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Recent developments have meant that this valuable technique can now be applied as an essential diagnostic tool in prenatal, haematological and pathological chromosomal analysis. Target DNA, after fixation and denaturation, is available for annealing to a similarly denatured, fluorescently labelled DNA probe, which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed and the DNA is counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

Probe Information

Cytocell's Chromoprobe Multiprobe® OctoChrome™ combines the utility of an 8 square Multiprobe device and the whole chromosome painting probe (labelled in 3 different colours), to allow all 24 chromosomes to be identified on a single slide. The OctoChrome™ device allows the simultaneous analysis of the whole genome on one slide in one hybridisation.

Probe Specification

Each square of the OctoChrome™ device carries whole chromosome painting probes for three different chromosomes in three different colour fluorophores, red, green and blue (Texas Red, FITC and Aqua/DEAC spectra respectively) which are visible simultaneously with a DAPI/FITC/Texas Red triple filter or specific single filters.

The arrangement of chromosome combinations on the OctoChrome™ device has been designed to facilitate the identification of the non-random chromosome rearrangements found in the most common leukaemias (Figure 1).

The OctoChrome™ is intended for FISH on metaphase chromosomes from fixed cultured peripheral blood cells.


Figure 1. Location of the probes on the Chromoprobe Multiprobe® OctoChrome™.

Materials Provided

Each kit contains the following reagents, which are sufficient for either 2 (PMP 802), 5 (PMP 804) or 10 (PMP 803) patient samples:

- 2, 5 or 10 Chromoprobe Multiprobe® OctoChrome™ devices coated with directly labelled probes:
Amount of red painting probe: 12-16ng per square
Amount of green painting probe: 36-51 ng per square
Amount of blue painting probe: 112-122ng per square
- 4, 7 or 12 glass slides printed with a special template
- 500µl Hybridisation Solution (Formamide, Dextran Sulphate, SSC)
- 1 Cytocell Slide Surface Thermometer
- 1 Cytocell Chromoprobe Multiprobe® Hybridisation Chamber

Warnings and Precautions

- For *in vitro* diagnostic use. For professional use only.
- Wear gloves when handling Hybridisation solution, DNA probes and DAPI counterstain.
- Hybridisation Solution contains formamide, which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
- DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
- Dispose of all hazardous materials according to your institution's guidelines for hazardous waste disposal.
- Operators must be capable of visually distinguishing between red, blue and green.
- Failure to adhere to the protocol may affect the performance and lead to false positive/negative results.

Storage and Handling

The Chromoprobe Multiprobe® kit should be stored at 2-8°C until the expiry date indicated on the kit label. Do not freeze.

Equipment and Materials Necessary but not Supplied

- 500µl Counterstain (DAPI antifade (ES: 0.125µg/ml DAPI (4,6-diamidino-2-phenylindole)).
- Hotplate (with a solid plate and accurate temperature control up to 80°C).
- Variable volume micropipettes and tips range 1µl - 200µl.
- Water bath with accurate temperature control at 72°C.
- Microcentrifuge tubes (0.5ml).
- Fluorescence microscope (Please see Fluorescence Microscope Recommendation section).
- Plastic or glass coplin jars.
- Forceps.
- Fluorescence grade microscope lens immersion oil.
- Bench top centrifuge.
- Microscope slides.
- Timer.
- 37°C water bath without stirrer.

Fluorescence Microscope Recommendation

Use a 100-watt mercury lamp and plan apochromat objectives x63 or x100 for optimal visualisation. Use a triple bandpass filter DAPI/FITC/Texas Red for optimal visualisation of the green and red fluorophores and DAPI simultaneously. The aqua fluorophore has specificity to the Aqua and DEAC spectrum (single bandpass Aqua or DEAC filter is required).

Check the fluorescence microscope before use to ensure it is operating correctly. Use immersion oil that is suitable for fluorescence microscopy and formulated for low autofluorescence. Avoid mixing DAPI Antifade with microscope immersion oil as this will obscure signals. Follow manufacturers' recommendations in regards to the life of the lamp and the age of the filters.

Sample Preparation

The Chromoprobe Multiprobe® OctoChrome™ is designed for use on cultured peripheral blood cells fixed in Carnoy's fixative that should be prepared according to the laboratory or institution guidelines.

Prepare air dried samples on Cytocell Chromoprobe Multiprobe® template slides according to Cytocell's protocol below. Baking or otherwise ageing slides is not recommended as it may reduce signal fluorescence.

Chromoprobe Multiprobe® Protocol


Please note: The probes used on the Chromoprobe Multiprobe® device are directly labelled with fluorophores, which are light sensitive. Ensure that exposure of the probes to laboratory lights is limited at all times (it is not necessary to work in the dark).

1. Slide preparation

- Clean a template slide. Soak the template slide for 2 minutes in 100% methanol and polish dry with a clean soft tissue.
- Establish the correct mitotic index. It is important that the intended sample has a sufficiently high mitotic index to allow detection of chromosome abnormalities. To check the density of the sample, using a micropipette (e.g. a Gilson P10 or P20) pipette 4µl of the cell suspension onto one of the areas of the spare template slide and allow to air dry. The small volume used means that you usually have to gently touch the slide with the pipette tip to transfer the suspension. Examine by phase contrast microscopy. If the cell density is too high, dilute the suspension with fresh fixative. If the mitotic index is too low, spin down the fixed cell suspension at 160xg for 10 minutes. Note the volume of supernatant, remove, and re-suspend the cell pellet in a smaller volume of fresh fixative. If cell sample density has been altered, spot 4µl of the concentrated sample onto another square of your test slide and re-examine by phase contrast microscopy.

Please Note: 50µl is the minimum volume required for the protocol.

- c) Quality control of samples. Samples should be examined for cytoplasm since this will interfere with the *in situ* protocol. If the chromosomes appear to be enclosed by a granular material when examined under phase contrast microscopy, then this will compromise results. One method for reducing cytoplasm is to spot 4µl of your sample onto the template slide and watch the fixative as it spreads out: in the normal situation, the fixative will spread to maximum, recede and then evaporate. To clean up any cytoplasm we have found that effective results are achieved if a fresh drop of fixative is allowed to fall onto the spot at the point when the spreading fixative has reached its maximum. Allow the drop of fixative to evaporate and re-examine the spot.
- d) Spotting of the slide. Pipette 4µl of cell suspension onto all 8 areas of the template slide in a sequence of alternating squares as shown below. This will prevent the cell spreads from interfering with each other.

- e) Once the first group of drops has air-dried, spot the remaining squares with 4µl drops in the same manner. After the slide has dried, examination of the slide under phase contrast will reveal whether any squares have been missed. If spots have been missed, or squares have too few cells, simply spot those squares again: it is not necessary to re-spot a new slide. If upon examination of the slide, a square has insufficient cells/metaphases, further drop(s) of suspension can be added to increase the cell density.

Please note: If the metaphases appear overspread, then clean a new template slide thoroughly in methanol and re-spot allowing every spot to dry before proceeding to the next.

2. Preparation of the Chromoprobe Multiprobe® device and template slide

- Ensure that the Chromoprobe Multiprobe® Hybridisation Chamber is in the 37°C water bath and allow it to equilibrate to 37°C (+/- 1°C). This may take up to an hour if the water bath has been switched on from cold.
- Mix the hybridisation solution by repeated pipetting and pre-warm a 25µl aliquot per device to 37°C. Also pre-warm each device by placing it on a 37°C hotplate label side down. Do not touch the raised surfaces of the device.
- Immerse template slides containing fixed samples in 2xSSC for 2 minutes at room temperature (RT) without agitation.
- Whilst the device is still at 37°C, dehydrate template slides containing fixed samples through an ethanol series (2 minutes each in 70%, 85% and 100%) at RT, air dry and place at 37°C hotplate to warm up.
- Add 2µl of pre-warmed hybridisation solution to each of the eight areas on the pre-warmed device using a P10 micropipette while it remains on a 37°C hotplate.

3. Positioning of template slide over the device

- Carefully invert the template slide over the device such that the number 1, which is now upside down, is located over the top right hand area of the device. To help locate square 1 (Chromosomes 3, 15 & 17), its position on the device has been marked with an orange dot).
- Make sure that the template slide is carefully aligned with the matching areas on the device. Carefully lower the slide over the device so that the drops of hybridisation solution make contact with the slide. Apply gentle, even pressure to ensure that the hybridisation solution is spread to the edges of each of the raised areas on the device.
- Lift the slide carefully holding the frosted end of the glass slide and invert so that the template slide is underneath the device. Make sure the device does not smear across the template slide as this could cause cross-contamination of the probes.
- Place at 37°C (+/- 1°C) (hotplate or incubator) for 10 minutes.

4. Instructions for use of the Cytocell Slide Surface Thermometer

- The temperature of the 75°C hotplate should be checked with the Cytocell Slide Surface Thermometer before proceeding to denaturation.
- To use the thermometer properly, place it onto the surface of the hotplate and wait until the different segments stop changing colour. The actual temperature is indicated by a deep aqua colour.

Please note:

- When the segments appear granular and the colours no longer appear uniform and regular, the thermometer should be discarded as it is exhausted. The life span of each thermometer should, however, easily be sufficient for a ten-device kit.
- This thermometer is a liquid crystal device and although reusable, it must be treated with care to ensure a reasonable life span. The thermometer must only be used to check the temperature of a hotplate; it must not be used to monitor the hotplate performance over time.

5. Denaturation

Please note: A PCR thermal cycler-heating block is NOT suitable for use in place of solid bed hotplate for this procedure.

- Transfer the slide/device sandwich to the hotplate taking particular care to hold it level. Ensure the sample slide is in good contact with the hotplate.
- Denature on the hotplate at 75°C (+/- 1°C) for 5 minutes.

6. Hybridisation

Place the slide/device sandwich in the pre-warmed Chromoprobe Multiprobe® Hybridisation Chamber, replace the lid and float the chamber in the 37°C (+/- 1°C) water bath (non-stirring) overnight.

Please note:

- Do not seal the lid on the hybridisation chamber.
- Do not place a lid on the water bath.
- Do not hybridise in an incubator.
- Please ensure that the hybridisation chamber is completely dry (i.e. no water or damp tissue inside the chamber).

The humidity inside the chamber is vital for optimal hybridisation. The correct levels will be achieved following those steps.

7. Post-hybridisation stringent washes

Please note: Avoid processing more than two slides through the stringency washes at any one time.

- Remove the device carefully from the slide.
- Immerse the slide in 0.4xSSC (pH 7.0) at 72°C (+/- 1°C) for 2 minutes without agitation.
- Drain the slide and immerse it in 2xSSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds without agitation.

8. Mounting and visualisation of results

- Drain the slide and apply 20µl of DAPI antifade to each end of the slide.
- Cover with a coverslip (24x50mm or 24x60mm), remove any bubbles and allow the colour to develop in the dark for 10 minutes.
- View with a fluorescence microscope.

Please note: Certain types of microscope have slide holders, which make it difficult to view the extreme ends of the slide. If this occurs then simply turn the slide through 180°, which will help with the viewing of the slide.

Stability of Finished Slides

FISHED slides remain analysable for up to 1 month if stored in the dark at/or below RT.

Procedural Recommendations

- Baking or ageing of slides may reduce signal fluorescence
- Hybridisation conditions may be adversely affected by the use of reagents other than those provided or recommended by Cytocell Ltd.
- Use a calibrated thermometer for measuring temperatures of solutions, waterbaths and incubators as these temperatures are critical for optimum product performance.
- The wash concentrations, pH and temperatures are important as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.
- Incomplete denaturation can result in lack of signal and over denaturation can also result in non-specific binding.
- Over hybridisation can result in additional or unexpected signals.
- Users should optimise the protocol for their own samples prior to using the test for diagnostic purposes.
- Slides should be independently scored by two analysts.
- Suboptimal conditions may result in non-specific binding that may be misinterpreted as a probe signal.

Expected Results

- The acrocentric chromosome painting probes contain short arm material. This is shared between the D & G group chromosomes, so cross-hybridisation in these regions may be observed.
- The chromosome 16 whole chromosome painting probe may show faint cross-hybridisation to the heterochromatic regions of chromosome Y.
- The Y chromosome painting probe contains the pseudoautosomal regions common with the X chromosome. Consequently, cross-hybridisation may be observed in these regions.
- Chromosomes 1, 5 and 19 have centromeric DNA sequences in common. Consequently cross-hybridisations in these regions may be observed between these chromosomes.
- Chromosome 1 whole chromosome painting probe may show a brighter signal at 1p36.
- Whole chromosome painting probes do not cover the Heterochromatic blocks in chromosomes 1, 9 and 16.

Limitations

Reporting and interpretation of FISH results should be consistent with professional standards of practice and should take into consideration other clinical and diagnostic information. This kit is intended as an adjunct to other diagnostic laboratory tests and therapeutic action should not be initiated on the basis of the FISH result alone. Failure to adhere to the protocol may affect the performance and lead to false positive/negative results. This kit has not been validated for purposes outside of the intended use stated.

Additional Information

For additional product information please contact the CytoCell Technical Support Department.

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Patents and Trademarks

Chromoprobe, Cytocell and Chromoprobe Multiprobe are registered trademarks of Cytocell Ltd.

The Chromoprobe principle is covered by international patents WO9314223, EP0623177. The design of the Multiprobe is a registered design and is also covered by a Design Patent No. 420,745.

This product contains technology licensed from Life Technologies Corporation and is available for human diagnostics or life science research use only.

Cytocell Ltd.
 3-4 Technopark
 Newmarket Road
 Cambridge, CB5 8PB, UK.
 T: +44(0)1223 294048
 F: +44(0)1223 294986
 E: probes@cytozell.com
 W: www.cytozell.com