**Instructions For Use**

REF: LPH 074-S / LPH 074 / LPH 074-20

**IGH/FGFR3 Plus Translocation, Dual Fusion Probe**

PROFESSIONAL USE ONLY

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
Further information available at [www.cytocell.com](http://www.cytocell.com)

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei from fixed cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Recent developments have meant that this valuable technique can now be applied as an essential diagnostic tool in prenatal, haematological and pathological chromosomal analysis. Target DNA, after fixation and denaturation, is available for annealing to a similarly denatured, fluorescently labelled DNA probe, which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed and the DNA is counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

**Probe Information**

The t(4;14)(p16.3;q32.3) translocation is found in 15-20% of multiple myeloma patients and results in the known dysregulation of two genes, one on each derivative chromosome: MMSET (Multiple Myeloma SET domain containing protein, also known as WHSC1) on the der(4) and FGFR3 (fibroblast growth factor receptor 3) on the der(14)<sup>1</sup>. The majority of the breakpoints on chromosome 4 occur within the MMSET locus resulting in a fusion gene where the 5' exons of MMSET are replaced with the VDJ region of the heavy chain locus. This results in MMSET being brought under the influence of the enhancer E $\mu$  on the der(4)<sup>1</sup>, FGFR3 lies at least 50kb from the breakpoints after translocation to chromosome 14 and FGFR3 expression is upregulated via proximity to the strong enhancer elements (particularly the Ea enhancer) in the immunoglobulin heavy chain<sup>2</sup>. Other studies state that the breakpoints on chromosome 4 lie in a 113kb region between FGFR3 and MMSET exon 5<sup>3</sup>. The breakpoint on chromosome 14 is almost exclusively in or near the switch region of the IGH locus. This t(4;14) translocation is often cytogenetically cryptic and was poorly described before the advent of FISH techniques.

**Probe Specification**  
FGFR3, 4p16.3, Red  
IGH, 14q32.33, Green

The IGH/FGFR3 *Plus* product consists of probes, labelled in green, proximal to the Constant, and within the Variable segment of the IGH region and FGFR3 probes, labelled in red. The FGFR3 probe mix contains a 118kb probe telomeric to FGFR3, including the D4S2561E marker and a second probe covering the 126kb region centromeric to MMSET, including the D4S1182 marker.

**Materials Provided**

**Probe:** 50 $\mu$ l per vial (5 tests), 100 $\mu$ l per vial (10 tests) or 200 $\mu$ l per vial (20 tests)  
Amount of red FGFR3 probe: 16-20ng/test  
Amount of green IGH probe: 42-52.5ng/test  
The probes are provided premixed in hybridisation solution (Formamide; Dextran Sulphate; SSC) and are ready to use.

**Counterstain:** 150 $\mu$ l per vial (15 tests) or 500 $\mu$ l per vial (50 tests)

The counterstain is DAPI antifade (ES: 0.125 $\mu$ g/ml DAPI (4,6-diamidino-2-phenylindole)).

**Warnings and Precautions**

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide, which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

**Storage and Handling**

The Aquarius<sup>®</sup> kit should be stored between -25°C to -15 °C in a freezer until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

**Equipment Necessary but not Supplied**

1. Hotplate (with a solid plate and accurate temperature control up to 80°C).
2. Variable volume micropipettes and tips range 1 $\mu$ l - 200 $\mu$ l.
3. Water bath with accurate temperature control at 72°C.
4. Microcentrifuge tubes (0.5ml).
5. Fluorescence microscope (Please see Fluorescence Microscope Recommendation section).
6. Plastic or glass coplin jars.
7. Forceps.
8. Fluorescence grade microscope lens immersion oil.
9. Bench top centrifuge.
10. Microscope slides.
11. 24x24mm coverslips.
12. Timer.
13. 37°C incubator.
14. Rubber solution glue.

**Fluorescence Microscope Recommendation**

For optimal visualisation of the probe we recommend a 100-watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all fluorophores and DAPI simultaneously.

**Sample Preparation**

The kit is designed for use on cultured peripheral blood cells or cultured bone marrow cells fixed in Carnoy's fixative that should be prepared according to the laboratory or institution guidelines. Prepare air dried samples on microscope slides according to standard cytogenetic procedures.

**FISH Protocol**

(Note: Please ensure that exposure of the probe to laboratory lights is limited at all times).

**Slide preparation**

1. Spot the cell sample onto a glass microscope slide. Allow to dry.
2. Immerse the slide in 2xSSC for 2 minutes at room temperature (RT) without agitation.
3. Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 minutes at RT.
4. Allow to dry.

**Pre-Denaturation**

5. Remove the probe from the freezer and allow it to warm to RT.
6. Ensure that the probe solution is uniformly mixed with a pipette.
7. Remove 10 $\mu$ l of probe per test, and transfer it to a microcentrifuge tube. Quickly return the remaining probe to the freezer.
8. Place the probe and the sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 5 minutes.
9. Spot 10 $\mu$ l of probe mixture onto the cell sample and carefully apply a coverslip. Seal with rubber solution glue and allow the glue to dry completely.

**Denaturation**

10. Denature the sample and probe simultaneously by heating the slide on a hotplate at 75°C (+/- 1°C) for 2 minutes.

**Hybridisation**

11. Place the slide in a humid, lightproof container at 37°C (+/- 1°C) overnight.

**Post-Hybridisation Washes**

12. Remove the coverslip and all traces of glue carefully.
13. Immerse the slide in 0.4xSSC (pH 7.0) at 72°C (+/- 1°C) for 2 minutes without agitation.
14. Drain the slide and immerse it in 2xSSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds without agitation.

- Drain the slide and apply 10µl of DAPI antifade onto each sample.
- Cover with a coverslip, remove any bubbles and allow the colour to develop in the dark for 10 minutes.
- View with a fluorescence microscope.

#### Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at/ or below RT.

#### Procedural Recommendations

- Baking or ageing of slides is not recommended as it may reduce signal fluorescence.
- Hybridisation conditions may be adversely affected by the use of reagents other than those provided or recommended by CytoCell Ltd.
- The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators as these temperatures are critical for optimum product performance.
- The wash concentrations, pH and temperatures are important as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.
- Incomplete denaturation can result in lack of signal and over denaturation can also result in non-specific binding.

#### Expected Results

In a normal cell these probes should appear as discrete red and green spots, one for each homologue (resulting in a 2R, 2G conformation). In a t(4;14)(p16.3;q32.3) patient there should be two yellow fusion signals in addition to the red and green signals of the normal chromosomes 4 and 14 respectively (1R, 1G, 2Y).

#### Known Cross-Reactivity

The green IGH probe may show faint cross hybridisation to a region on chromosome 16.

#### Limitations

Reporting and interpretation of FISH results should be consistent with professional standards of practice and should take into consideration other clinical and diagnostic information. This kit is intended as an adjunct to other diagnostic laboratory tests and therapeutic action should not be initiated on the basis of the FISH result alone.

#### Additional Information

For additional product information please contact the CytoCell Technical Support Department.

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**Patents and Trademarks**
Aquarius and CytoCELL are registered trademarks of CytoCELL Ltd. This product contains technology licensed from Life Technologies Corporation and is available for human diagnostics or life science research use only.
**CytoCELL Ltd.**