**Instructions For Use**  
REF: LPS 049-S / LPS 049

**FOXO1 Breakpart Probe**

**PROFESSIONAL USE ONLY**

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Further information available at [www.cytocell.com](http://www.cytocell.com)

**Intended Use**

This product is intended to identify, via fluorescence *in situ* hybridisation (FISH), rearrangements involving the FOXO1 region on chromosome 13 at location 13q14.1 in patients with solid tissue neoplasia. It is intended to be used on formalin-fixed paraffin-embedded (FFPE) tissues. This product is for professional use only and is intended to be an adjunct to clinical cytogenetics.

**Background**

Fluorescence *in situ* hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei from fixed cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Recent developments have meant that this valuable technique can now be applied to the assessment of solid tumour biopsies as well, which can provide important information for prediction of tumour progression. Current methodologies, namely immunohistochemistry or blotting can provide information at the gene expression level but, when carried out on tissue sections, FISH can provide information at a gene level, *in situ*, at the precise site within the tumour. This can reveal cell-to-cell heterogeneity and enable the detection of small clones of genetically distinct cells.

**Probe Information**

Translocations involving the FOXO1 (*forkhead box O1*) gene at 13q14.1 and either the PAX3 (*paired box 3*) gene at 2q36.1 or the PAX7 (*paired box 7*) gene at 1p36.1 are seen frequently in cases of alveolar rhabdomyosarcoma<sup>1,2</sup>. Rhabdomyosarcoma is the most common soft-tissue sarcoma seen in children and younger adults with two major histological subtypes: alveolar rhabdomyosarcoma (ARMS) and embryonal rhabdomyosarcoma (ERMS)<sup>3</sup>. FOXO1 rearrangements are recognised recurrent abnormalities seen in ARMS, but not seen in ERMS<sup>1,2</sup>.


Approximately 55% of cases of ARMS are associated with a PAX3-FOXO1 rearrangement via a t(2;13)(q36.1;q14) translocation and 22% of cases of ARMS are associated with a PAX7-FOXO1 rearrangement via a t(1;13)(p36;q14) translocation<sup>4</sup>. These translocations lead to the fusion of transcription factor FOXO1 to the transcription factors PAX3 and PAX7 located at 2q36.1 and 1p36.13 respectively<sup>2</sup>.

Studies have shown that ARMS patients with PAX-FOXO1 rearrangements have an inferior outcome compared to ERMS patients, whereas ARMS patients without PAX-FOXO1 rearrangements show similar outcomes to ERMS<sup>2,5</sup>.

A subset of patients with ARMS may show fusion gene amplification. This is most commonly associated with the presence of PAX7-FOXO1 rearrangements and has been shown to be associated with significantly improved outcome over ARMS patients with PAX-FOXO1 rearrangements without fusion gene amplification<sup>6</sup>. This breakpart probe design allows the detection of FOXO1 rearrangements, irrespective of the partner gene involved.

**Probe Specification**

FOXO1, 13q14.1, Red  
FOXO1, 13q14.1, Green

The FOXO1 Breakpart probe consists of two probes (405kb and 183kb), labelled in green, situated proximal to the FOXO1 gene and covering markers D13S765

and SHGC-111293 and four probes (123kb, 142kb, 95kb and 253kb), labelled in red, situated distal to the FOXO1 gene and covering markers D13S638 and SHGC-16596.

**Materials Provided**

**Probe:** 50µl per vial (5 tests) or 100µl per vial (10 tests)

Amount of Red FOXO1 probe: 34-59ng/test

Amount of Green FOXO1 probe: 342.5-512.5ng/test

The probes are premixed in hybridisation solution (Formamide; Dextran Sulphate; SSC) and are ready to use.

**Counterstain:** 150µl per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125µg/ml DAPI (4,6-diamidino-2-phenylindole)).

**Warnings and Precautions**

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide, which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. Dispose of all hazardous materials according to your institution's guidelines for hazardous waste disposal.
6. Operators must be capable of visually distinguishing between red, blue and green.
7. Failure to adhere to the protocol may affect the performance and lead to false positive/negative results.

**Storage and Handling**

The Aquarius<sup>®</sup> kit should be stored between -25°C to -15 °C in a freezer until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark. Ensure that exposure of the probe and counterstain to laboratory lights is limited at all times.

**Equipment Necessary but not Supplied**

1. Hotplate (with a solid plate and accurate temperature control up to 80°C).
2. Variable volume micropipettes and tips range 1µl - 200µl.
3. Water bath with accurate temperature control at 72°C.
4. Microcentrifuge tubes (0.5ml).
5. Fluorescence microscope (Please see Fluorescence Microscope Recommendation section).
6. Plastic or glass coplin jars.
7. Forceps.
8. Fluorescence grade microscope lens immersion oil.
9. Bench-top centrifuge.
10. Microscope slides.
11. 24x24mm coverslips.
12. Timer.
13. 37°C incubator.
14. Rubber solution glue.
15. Aquarius<sup>®</sup> Tissue Pretreatment Kit (LPS 100).

**Fluorescence Microscope Recommendation**

Use a 100-watt mercury lamp and plan apochromat objectives x63 or x100 for optimal visualisation. Use a triple bandpass filter DAPI/FITC/Texas Red for optimal visualisation of the green and red fluorophores and DAPI simultaneously. Alternatively for viewing red and green fluorophores use dual bandpass filter FITC/Texas Red.

Check the fluorescence microscope before use to ensure it is operating correctly. Use immersion oil that is suitable for fluorescence microscopy and formulated for low autofluorescence. Follow manufacturers' recommendations in regards to the life of the lamp and the age of the filters.

**Sample Preparation**

The kit is designed for use on Formalin Fixed Paraffin Embedded (FFPE) tissue sections and should be prepared according to the laboratory or institution guidelines. Use 4µm - 6µm thick FFPE tissue sections for FISH.

**Tissue Sample Pretreatment**

Tissue sample pretreatment should be done according to the laboratory or institution guideline. Use the Aquarius<sup>®</sup> Tissue Pretreatment Kit (LPS 100) for optimal results.

**FISH Protocol**

(Note: Ensure that exposure of the probe and counterstain to laboratory lights is limited at all times)

**Pre-Denaturation**

1. Remove the probe from the freezer and allow it to warm to room temperature (RT). Briefly centrifuge tubes before use
2. Ensure that the probe solution is uniformly mixed with a pipette.
3. Remove 10µl of probe per test, and transfer it to a microcentrifuge tube. Quickly return the remaining probe to the freezer.
4. Place the probe and the sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 5 minutes.
5. Spot 10µl - 15µl of probe mixture onto the sample and carefully apply a coverslip. Seal with rubber solution glue and allow the glue to dry completely.

**Denaturation**

6. Denature the sample and probe simultaneously by heating the slide on a hotplate at 75°C (+/- 1°C) for 5 minutes.

**Hybridisation**

7. Place the slide in a humid, lightproof container at 37°C (+/- 1°C) overnight.

## Post-Hybridisation Washes

- Remove the coverslip and all traces of glue carefully.
- Immerse the slide in 0.4xSSC (pH 7.0) at 72°C (+/- 1°C) for 2 minutes without agitation.
- Drain the slide and immerse it in 2xSSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds without agitation.
- Drain the slide and apply 10µl - 15µl of DAPI antifade onto each sample.
- Cover with a coverslip, remove any bubbles and allow the colour to develop in the dark for 10 minutes.
- View with a fluorescence microscope.

## Comments

Hybridisation efficiency and tissue morphology are usually negatively correlated. Aggressive pretreatment procedures improving hybridisation efficiency (e.g. an extended enzyme digestion time) tend to destroy cell structure and tissue morphology. However, mild pretreatment saving tissue structures may not be sufficient for probe penetration and successful FISH results.

The optimal length of heat pretreatment and enzyme digestion time will depend on the age of the block, the tissue composition, and quality of tissue fixation. Decrease enzyme digestion for core biopsies and any sections that contain few tumour cells or have large areas of necrosis. Handle these samples particularly carefully to avoid over-digestion.

## Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark below 4°C.

## Procedural Recommendations

- Hybridisation conditions may be adversely affected by the use of reagents other than those provided or recommended by CytoCell Ltd.
- Use a calibrated thermometer for measuring temperatures of solutions, waterbaths and incubators as these temperatures are critical for optimum product performance.
- The wash concentrations, pH and temperatures are important as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.
- Incomplete denaturation can result in lack of signal and over denaturation can also result in non-specific binding.
- Over hybridisation can result in additional or unexpected signals.
- Users should optimise the protocol for their own samples prior to using the test for diagnostic purposes.

## Expected Results

In a normal cell, two red/green fusion signals are expected (2F). In a cell with a balanced FOXO1 rearrangement, the expected signal pattern will be 1 fusion, 1 red and 1 green (1R, 1G, 1F). In a cell with amplification of proximal FOXO1 the expected signal pattern will be 1 fusion, 1 red and amplified green. Other signal patterns are possible in aneuploid/unbalanced specimens.

## Known Cross-Reactivity

No known cross-reactivity.

## Limitations

Reporting and interpretation of FISH results should be consistent with professional standards of practice and should take into consideration other clinical and diagnostic information. This kit is intended as an adjunct to other diagnostic laboratory tests and therapeutic action should not be initiated on the basis of the FISH result alone.

Failure to adhere to the protocol may affect the performance and lead to false positive/negative results.

This kit has not been validated for purposes outside of the intended use stated.

## Additional Information

For additional product information please contact the CytoCell Technical Support Department.

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#### Patents and Trademarks
Aquarius and CytoCELL are registered trademarks of CytoCELL Ltd.
**CytoCELL Ltd.**