

PROBLEM	INSUFFICIENT FIXATION ON BLOOD SMEARS	HIGH LEVELS OF GREEN BACKGROUND	BRIGHT, GREEN SIGNALS ON CELL SEPARATED SAMPLES	SMALL, WEAK SIGNALS	BACTERIAL CONTAMINATION
<p>LIKELY CAUSES</p>	<ul style="list-style-type: none"> <li>High levels of red background caused by presence of erythrocytes and cellular debris on the slide</li> <li>The fixation time is critical and depends on the size and erythrocyte density of the smear</li> </ul>	<ul style="list-style-type: none"> <li>Positively-charged slides are designed for tissue sections to ensure good adhesions, however these can produce higher levels of green background on other cell samples</li> </ul>	<ul style="list-style-type: none"> <li>CD138+ separated samples can produce a high volume of cytoplasm and debris, which can be seen as bright green fluorescence</li> </ul>	<ul style="list-style-type: none"> <li>Caused by baking or ageing slides</li> <li>Could also be a result of high stringency in the post-hybridisation wash</li> </ul>	<ul style="list-style-type: none"> <li>Bacteria may fluoresce very brightly under both the FITC (Fluorescein isothiocyanate) and Texas Red® filters, depending on the type of bacteria</li> <li>The level of contamination can vary depending on the length of culture i.e 24 or 72 hours</li> </ul>
<p>SOLUTION</p>	<ul style="list-style-type: none"> <li>✓ 10-30 minutes fixation time with 3:1 methanol/ glacial acetic acid is recommended</li> <li>✓ Replace with fresh fixative (3:1 methanol/ glacial acetic acid fix) after 10 minutes</li> <li>✓ Hypotonic solutions, such as potassium chloride can be used with smear slides to improve FISH results</li> </ul>	<ul style="list-style-type: none"> <li>✓ Review the type of slides used, try using non-adhesive slides</li> </ul>	<ul style="list-style-type: none"> <li>✓ Use &lt;3µl cell suspension and treat with a hypotonic solution, such as potassium chloride, followed by an enzymatic digestion step using a protease</li> </ul>	<ul style="list-style-type: none"> <li>✓ Avoid baking (&gt;45°C) or ageing slides in heated SSC solutions. This reduces signal intensity and can be time consuming</li> <li>✓ Re-make 0.4xSSC ensuring that the pH=7 and the wash temperature is 72+/-1°C</li> </ul>	<ul style="list-style-type: none"> <li>✓ Always use fresh cell samples and reset the culture if required</li> <li>✓ Ensure aseptic techniques are correctly followed when culturing cells to minimise risk of contamination</li> <li>✓ Regularly refresh your solutions, including reagent bottles and coplin jars</li> </ul>

## INTERPRETATION OF SIGNALS

**Breakapart:**

Normal Signal Pattern (2F)

Abnormal Signal Pattern (1R, 1G, 1F)

In some breakapart probes, the red and green signals may not be tightly co-localised in all cells, especially when chromatin is extended in interphase. It is therefore important to determine the cut-off levels for a normal result<sup>1</sup>.

**Translocation:**

Normal Signal Pattern (2R, 2G)

Abnormal Signal Pattern (1R, 1G, 2F)

The dual-fusion translocation probes span the breakpoints of interest. The abnormal pattern is rarely seen as an artefact of chromosome overlap in interphase<sup>1</sup>. In an abnormal cell, this produces two fusion signals (co-localised red and green signals).

**Scoring Recommendations<sup>1</sup>:**

- Analyse only intact nuclei, not overlapped or crowded nuclei or nuclei covered by cytoplasmic debris or high degree of auto-fluorescence
- Signal intensity may vary, even within a single nucleus. In such cases, use single filters and/or adjust the focal plane.

1. Analysis and scoring guidelines recommended by The AGT Cytogenetics Laboratory Manual (fourth edition).



For more troubleshooting tips, visit our FISH 'n' Tips page:  
[www.cytocell.com/support](http://www.cytocell.com/support)  
[contact@ogt.com](mailto:contact@ogt.com)