A Sysmex Group Company

Instructions For Use  
REF: LPS 018-S / LPS 018

## FGFR1 Breakpart/Amplification Probe

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Further information available at [www.cytocell.com](http://www.cytocell.com)


Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei from fixed cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Recent developments have meant that this valuable technique can now be applied to the assessment of solid tumour biopsies as well, which can provide important information for prediction of tumour progression. Current methodologies, namely immunohistochemistry or Southern blotting can provide information at the gene expression level but, when carried out on tissue sections (either cryostat or paraffin embedded), FISH can provide information at a gene level, *in situ*, at the precise site within the tumour. This can reveal cell-to-cell heterogeneity and enable the detection of small clones of genetically distinct cells.

### Probe Information

Fibroblast Growth Factor Receptor 1 (FGFR1), located on chromosome 8, was the first gene to be shown to be amplified in some breast cancers<sup>1</sup> and is also translocated in some patients with 8p11 myeloproliferative syndrome<sup>2</sup>. Amplification of the gene is thought to be the causative factor in around 10% of breast tumours<sup>3</sup> and has been shown to provide a poor prognosis to patients as over-expression of the gene product can lead to early relapse<sup>4</sup>. In common with other genetic defects in lymphoproliferative disorders, like the BCR/ABL translocation in CML, the FGFR1 gene contains sequences that could code for a putative tyrosine kinase protein<sup>5,6</sup>, which has been implicated in myeloproliferative syndrome (also known as EMS or Stem cell leukaemia/lymphoma syndrome, SCLL), marking out the importance of this gene in both solid tumours and cancers of the blood. There are at least eight partner genes involved with FGFR1 translocations, including ZNF198 (the most common<sup>7</sup>), FOP, CEP110, BCR, HERV-K, FGFR10P2, TIF1 and MYO18A<sup>8</sup>. Cytocell has developed a three-colour combined breakpart and amplification FISH probe for FGFR1 that can be used in either bone marrow samples from myeloproliferative syndrome (EMS) patients or tissue sections from patients where gene amplification may be suspected. In EMS, the breakpart strategy will show the split of one of the two fusion signals along with a blue centromere to enumerate chromosome 8. In patients that are being tested for amplifications of FGFR1, the non-translocated FGFR1 FISH probe will appear as a fusion signal and this will appear amplified. The chromosome 8 centromere probe will also act as a control probe in these cases.

### Probe Specification

FGFR1, 8p11.23-p11.22, Red  
FGFR1, 8p11.23-p11.22, Green  
D8Z2, 8p11.1-q11.1, Blue

The FGFR1 Breakpart/Amplification probe consists of a green 272kb probe and a red 267kb probe, which are positioned on each side of the FGFR1 gene. The 8-centromere probe in blue acts as a control for chromosome 8.

### Materials Provided

**Probe:** 50µl per vial (5 tests) or 100µl per vial (10 tests)  
Amount of Red FGFR1 probe: 80-100ng/test  
Amount of Green FGFR1 probe: 240-300ng/test  
Amount of Blue 8 cen probe: 84-105ng/test  
The probes are provided premixed in hybridisation solution (Formamide; Dextran Sulphate; SSC) and are ready to use.

### Counterstain: 150µl per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125µg/ml DAPI (4,6-diamidino-2-phenylindole)).

### Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide, which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

### Storage and Handling

The Aquarius<sup>®</sup> Kit should be stored between -25°C to -15 °C in a freezer until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

### Equipment Necessary but not Supplied

1. Hotplate (with a solid plate and accurate temperature control up to 80°C).
2. Variable volume micropipettes and tips range 1µl - 200µl.
3. Water bath with accurate temperature control at 72°C.
4. Microcentrifuge tubes (0.5ml).
5. Fluorescence microscope (Please see Fluorescence Microscope Recommendation section).
6. Plastic or glass coplin jars.
7. Forceps.
8. Fluorescence grade microscope lens immersion oil.
9. Bench-top centrifuge.
10. Microscope slides.
11. 24x24mm coverslips.
12. Timer.
13. 37°C incubator.
14. Rubber solution glue.
15. Aquarius<sup>®</sup> Tissue Pretreatment Kit (LPS 100).

### Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100-watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all fluorophores and DAPI simultaneously. Alternatively for viewing red and green fluorophores use dual bandpass filter FITC/Texas Red. The blue fluorophore has specificity to the Aqua and DEAC spectrum (single bandpass Aqua or DEAC filter is required).

### Sample Preparation

This Kit has been designed for use on:

- Formalin Fixed Paraffin Embedded (FFPE) Tissue sections or Tissue Microarrays (TMA), 4µm - 6µm thick tissue sections should be used.
- Peripheral blood samples or cultured bone marrow cells fixed in Carnoy's fixative air dried on a microscope slides according to standard cytogenetic procedures.

All samples should be prepared according to the laboratory or institution guidelines.

### FISH Protocol

(Note: Please ensure that exposure of the probe to laboratory lights is limited at all times)

### FFPE Procedure

#### Tissue Sample Pretreatment

Tissue sample pretreatment should be done according to the laboratory or institution guideline. For optimal results use the Aquarius<sup>®</sup> Tissue Pretreatment Kit (LPS 100).

### Pre-Denaturation

1. Remove the probe from the freezer and allow it to warm to room temperature (RT)
2. Ensure that the probe solution is uniformly mixed with a pipette.
3. Remove 10µl - 15µl (depending on the size of the tissue) of probe per test, and transfer it to a microcentrifuge tube. Quickly return the remaining probe to the freezer.
4. Place the probe and the sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 5 minutes.
5. Spot 10µl - 15µl of probe mixture onto the sample and carefully apply a coverslip. Seal with rubber solution glue and allow the glue to dry completely.

### Denaturation

6. Denature the sample and probe simultaneously by heating the slide on a hotplate at 75°C (+/- 1°C) for 5 minutes.

## Hybridisation

- Place the slide in a humid, lightproof container at 37°C (+/- 1°C) overnight.

## Post-Hybridisation Washes

- Remove the coverslip and all traces of glue carefully.
- Immerse the slide in 0.4xSSC (pH 7.0) at 72°C (+/- 1°C) for 2 minutes without agitation.
- Drain the slide and immerse it in 2xSSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds without agitation.
- Drain the slide and apply 10µl - 15µl of DAPI antifade onto each sample.
- Cover with a coverslip, remove any bubbles and allow the colour to develop in the dark for 10 minutes.
- View with a fluorescence microscope.

## Peripheral blood or cultured bone marrow Procedure

### Slide preparation

- Spot the cell sample onto a glass microscope slide. Allow to dry.
- Immerse the slide in 2xSSC for 2 minutes at room temperature (RT) without agitation.
- Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 minutes at RT.
- Allow to dry.

### Pre-Denaturation

- Remove the probe from the freezer and allow it to warm to RT.
- Ensure that the probe solution is uniformly mixed with a pipette.
- Remove 10µl of probe per test, and transfer it to a microcentrifuge tube. Quickly return the remaining probe to the freezer.
- Place the probe and the sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 5 minutes.
- Spot 10µl of probe mixture onto the cell sample and carefully apply a coverslip. Seal with rubber solution glue and allow the glue to dry completely.

### Denaturation

- Denature the sample and probe simultaneously by heating the slide on a hotplate at 75°C (+/- 1°C) for 2 minutes.

### Hybridisation

- Place the slide in a humid, lightproof container at 37°C (+/- 1°C) overnight.

### Post-Hybridisation Washes

- Remove the coverslip and all traces of glue carefully.
- Immerse the slide in 0.4xSSC (pH 7.0) at 72°C (+/- 1°C) for 2 minutes without agitation.
- Drain the slide and immerse it in 2xSSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds without agitation.
- Drain the slide and apply 10µl - 15µl of DAPI antifade onto each sample.
- Cover with a coverslip, remove any bubbles and allow the colour to develop in the dark for 10 minutes.
- View with a fluorescence microscope.

### Comments

Hybridisation efficiency and tissue morphology are usually negatively correlated. Aggressive pretreatment procedures improving hybridisation efficiency (e.g. an extended enzyme digestion time) tend to destroy cell structure and tissue morphology. However, mild pretreatment saving tissue structures may not be sufficient for probe penetration and successful FISH results.

The optimal length of heat pretreatment and enzyme digestion time will depend on the age of the block, the tissue composition, and quality of tissue fixation. Enzyme digestion should be decreased for core biopsies and any sections that contain few tumour cells or have large areas of necrosis. These samples need to be handled particularly carefully to avoid over-digestion.

### Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark below 4°C.

### Procedural Recommendations

- Baking or ageing of slides containing peripheral blood samples or bone marrow is not recommended as it may reduce signal fluorescence.
- Hybridisation conditions may be adversely affected by the use of reagents other than those provided or recommended by Cytocell Ltd.
- The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators, as these temperatures are critical for optimum product performance.
- The wash concentrations, pH and temperatures are important as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.
- Incomplete denaturation can result in lack of signal and over denaturation can also result in non-specific binding.

### Expected Results

Cytocell has developed a three-colour combined breakapart and amplification FISH probe for FGFR1, which can be used in either bone marrow samples from EMS patients, or tissue section in possibly amplified patients.

As it is the involvement of the FGFR1 tyrosine kinase domain that is important in EMS, detection of the partner gene involved is less important than establishing that there has been a translocation, which has disrupted the FGFR1 gene. In these cases, the breakapart strategy will show the conventional split of one of the two fusion signals into a red and a green signal, along with a blue 8 centromere to enumerate chromosome 8. Although not strictly relevant in these cases, this additional probe may act as an extra enumerator in some cases.

In patients that are being tested for amplifications of FGFR1, the breakapart strategy is not relevant, but the un-translocated FGFR1 FISH signal will appear as a fusion probe and it is this fusion signal that will appear amplified. In order to

properly enumerate this probe in tissue sections, the chromosome 8 centromere probe, included in a contrasting colour (blue), is useful.

In a normal cell due to the close juxtaposition of the probes, 2 yellow fusion and 2 blue control signals will be observed (2Y, 2B). In the event of a 8p11 translocation one split and one fused signal will be observed resulting in 1 yellow, 1 red, 1 green, 2 blue conformation (1Y, 1R, 1G, 2B). In the event of an amplification of the FGFR1 gene 3 or more yellow fusion and 2 blue control signals will be observed (3+Y, 2B).

### Limitations

Reporting and interpretation of FISH results should be consistent with professional standards of practice and should take into consideration other clinical and diagnostic information. This kit is intended as an adjunct to other diagnostic laboratory tests and therapeutic action should not be initiated on the basis of the FISH result alone.

### Additional Information

For additional product information please contact the Cytocell Technical Support Department.

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#### Patents and Trademarks
Aquarius and CytoCELL are registered trademarks of CytoCELL Ltd.
**CytoCELL Ltd.**